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
## BTS 310<sup>®</sup>

REVISÃO: 15 (11/2012)

## ÍNDICE

REVISÃO: 15 (11/2012)

KIT	PÁGINA	KIT	PÁGINA
ALBUMINA .....	01	FOSFATASES .....	23
ALDOLASE .....	02	FOSFATO .....	24
ALT/TGP .....	03	FOSFATO UV .....	25
AMILASE .....	04	GAMAGT .....	26
AMILASE CNPG .....	05	GAMA GT CINÉTICA .....	27
AST/TGO .....	06	GLICOSE ENZ. LÍQUIDA .....	28
BILIRRUBINA .....	07	GLUCOX 500 .....	29
CÁLCIO .....	08	HEMOGLOBINA .....	30
CÁLCIO ARSENAZO .....	09	HEMOGLOBINA GLICADA .....	31
CKMB (LANÇAMENTO) .....	10	MAGNÉSIO .....	32
CK NAC .....	11	MICROPROTE pirogalol .....	33
CLORETO COLORIMÉTRICO .....	12	MUCOPROTEÍNAS .....	34
COLESTEROL ENZ. LÍQUIDO .....	13	POTÁSSIO .....	35
COLESTEROL HDL .....	14	PROTEÍNAS TOTAIS .....	36
COLESTEROL 250.....	15	TRANSAMINASES .....	37
COLINESTERASE .....	16	TRIGLICÉRIDES 120 .....	38
CPK.....	17	TRIGLICÉRIDES ENZ. LÍQUIDO .....	39
CREATININA CINÉTICA .....	18	URATO 160 .....	40
DHL.....	19	URATO ENZIMÁTICO LÍQUIDO .....	41
DHL UV .....	20	UREIA 500 .....	42
FERRO .....	21	UREIA UV .....	43
FOSFATASE ALCALINA CINÉTICA .....	22		

## ALBUMINA

### END POINT WITH STANDARD

UNITS:	g/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	670
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	4.000

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## ALT/TGP

### DELTA KINETICS WITH FACTOR

UNIT:	UI/L
READING FILTER:	340nm
FACTOR:	1746
STABILIZ. TIME:	1
INCUBATION TIME:	60
INTERVAL TIME:	20
INTERVALS NUMBER:	4
SIPPING VOLUME:	400
REACTION SLOPE:	Decrease
TEMPERATURE:	37

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## AMILASE CNPG

### DELTA KINETICS WITH FACTOR

UNIT:	UI/L
READING FILTER:	405nm
FACTOR:	3953
STABILIZ. TIME:	1
INCUBATION TIME:	30
INTERVAL TIME:	60
INTERVALS NUMBER:	2
SIPPING VOLUME:	500
REACTION SLOPE:	Increase
TEMPERATURE:	37

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#### . Bilirrubina Total

DIFERENTIAL MODE WITH STANDARD

UNITS:	mg/dL
READING FILTER:	546
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
STANDARD:	*

#### . Bilirrubina direta

DIFERENTIAL MODE WITH STANDARD

UNITS:	mg/dL
READING FILTER:	546
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
STANDARD:	*

### . TÉCNICA MICRO

#### . Bilirrubina Total

DIFERENTIAL MODE WITH STANDARD

UNITS:	mg/dL
READING FILTER:	546
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
STANDARD:	*

#### . Bilirrubina direta

DIFERENTIAL MODE WITH STANDARD

UNITS:	mg/dL
READING FILTER:	546
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
STANDARD:	*

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## CÁLCIO

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	578
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	*
REACTION SLOPE:	INCREASING
STANDARD:	10.00

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## CÁLCIO ARSENAZO

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	670
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	500
REACTION SLOPE:	INCREASING
STANDARD:	10

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## CK MB

### FIXED TIME WITH FACTOR

UNIT:	U./L
READING FILTER:	340nm
FACTOR:	1350
STABILIZ. TIME:	1
INCUBATION TIME:	180
INTERVAL TIME:	300
SIPPING VOLUME:	500
REACTION SLOPE:	Increase
TEMPERATURE:	37

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## CK-NAC

### DELTA KINETICS WITH FACTOR

UNIT:	U/L
READING FILTER:	340nm
FACTOR:	8199
STABILIZ. TIME:	1
INCUBATION TIME:	120
INTERVAL TIME:	60
INTERVALS NUMBER:	3
SIPPING VOLUME:	400
REACTION SLOPE:	Increase
TEMPERATURE:	37

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## CLORETOS COLORIMÉTRICO

### END POINT WITH STANDARD

UNITS:	mEq/L
READING MODE:	MONOCHROMATIC
READING FILTER:	505
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	100.0

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## COLESTEROL ENZ. LÍQUIDO

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	505
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	200.0

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## COLESTEROL HDL

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	505
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	100.0

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## CREATININA CINÉTICA

### FIXED TIME WITH STANDARD

UNIT:	mg/dL
READING FILTER:	505nm
STABILIZ. TIME:	1
INCUBATION TIME:	30
INTERVAL TIME:	60
SIPPING VOLUME:	500
REACTION SLOPE:	Increase
STANDARD:	5
TEMPERATURE:	37

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## DHL - UV

### DELTA KINETICS WITH FACTOR

UNIT:	UI/L
READING FILTER:	340nm
FACTOR:	8109
STABILIZ. TIME:	1
INCUBATION TIME:	30
INTERVAL TIME:	60
INTERVALS NUMBER:	3
SIPPING VOLUME:	400
REACTION SLOPE:	Decrease
TEMPERATURE:	37

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## FOSFATASE ALCALINA CINÉTICA

### DELTA KINETICS WITH FACTOR

UNIT:	U/L
READING FILTER:	405nm
FACTOR:	2764
STABILIZ. TIME:	1
INCUBATION TIME:	60
INTERVAL TIME:	60
INTERVALS NUMBER:	3
SIPPING VOLUME:	400
REACTION SLOPE:	Increase
TEMPERATURE:	37

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## FOSFATO

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	670
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	5.000

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## FOSFATO - UV

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	340
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	5.000

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## GAMA GT CINÉTICA

### DELTA KINETICS WITH FACTOR

UNIT:	U/L
READING FILTER:	405nm
FACTOR:	1111
STABILIZ. TIME:	1
INCUBATION TIME:	20
INTERVAL TIME:	60
INTERVALS NUMBER:	3
SIPPING VOLUME:	400
REACTION SLOPE:	Increase
TEMPERATURE:	37

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## MAGNÉSIO

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	505
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	2.000

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## MICROPROTE pirogalol

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	578
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	500
REACTION SLOPE:	INCREASING
STANDARD:	100

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## MUCOPROTEÍNAS

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	670
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	6.000

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## POTÁSSIO

### END POINT WITH STANDARD

UNITS:	mEq/L
READING MODE:	MONOCHROMATIC
READING FILTER:	578
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	5.000

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## PROTEÍNAS TOTAIS

### END POINT WITH STANDARD

UNITS:	g/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	546
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	4.000

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.TGO/AST

MULTISTANDARD SPLINE

UNITS: F.R./mL  
READING MODE: MONOCHROMATIC  
READING FILTER: 505  
STABILIZ. TIME: 3  
TEMPERATURE: 37  
SIPPING VOLUME: 900  
REACTION SLOPE: INCREASING  
REPLICATES STD.: 1  
REPLICATES SAMP.: 1

AXIS

ABS: LIN  
CONC: LIN  
STANDARD 1: 0.000  
STANDARD 2: 24.00  
STANDARD 3: 61.00  
STANDARD 4: 114.0  
STANDARD 5: 190.0  
DILUTION FACTOR: 1.000

.TGP/ALT

MULTISTANDARD SPLINE

UNITS: F.R./mL  
READING MODE: MONOCHROMATIC  
READING FILTER: 505  
STABILIZ. TIME: 3  
TEMPERATURE: 37  
SIPPING VOLUME: 900  
REACTION SLOPE: INCREASING  
REPLICATES STD.: 1  
REPLICATES SAMP.: 1

AXIS

ABS: LIN  
CONC: LIN  
STANDARD 1: 0.000  
STANDARD 2: 28.00  
STANDARD 3: 57.00  
STANDARD 4: 97.00  
STANDARD 5: 150.0  
DILUTION FACTOR: 1.000

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.TGO/AST

MULTISTANDARD SPLINE

UNITS: U.I./L  
READING MODE: MONOCHROMATIC  
READING FILTER: 505  
STABILIZ. TIME: 3  
TEMPERATURE: 37

SIPPING VOLUME: 900  
REACTION SLOPE: INCREASING  
REPLICATES STD.: 1  
REPLICATES SAMP.: 1  
AXIS  
ABS: LIN  
CONC: LIN  
STANDARD 1: 0.000  
STANDARD 2: 11.50  
STANDARD 3: 29.40  
STANDARD 4: 54.94  
STANDARD 5: 91.58  
DILUTION FACTOR: 1.000

.TGP/ALT

MULTISTANDARD SPLINE

UNITS: U.I./L  
READING MODE: MONOCHROMATIC  
READING FILTER: 505  
STABILIZ. TIME: 3  
TEMPERATURE: 37  
SIPPING VOLUME: 900  
REACTION SLOPE: INCREASING  
REPLICATES STD.: 1  
REPLICATES SAMP.: 1

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## URATO 160

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	505
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	7.000

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## UREIA 500

### END POINT WITH STANDARD

UNITS:	mg/dL
READING MODE:	MONOCHROMATIC
READING FILTER:	670
STABILIZ. TIME:	3
TEMPERATURE:	37
SIPPING VOLUME:	900
REACTION SLOPE:	INCREASING
STANDARD:	80.00

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## UREIA-UV

### FIXED TIME WITH STANDARD

UNIT:	mg/dL
READING FILTER:	340nm
STABILIZ. TIME:	1
INCUBATION TIME:	30
INTERVAL TIME:	60
SIPPING VOLUME:	400
REACTION SLOPE:	Decrease
STANDARD:	80
TEMPERATURE:	37

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