AUTOHUMALYSER 900 S[®]

LANÇAMENTOS!

**CÁLCIO ARSENAZO
MICROPROTE pirogalol**

REVISÃO: 04 (02/2009)

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TEMPERATURA: 37°C

TESTNAME	ALBUMINA
DECIMAL	0
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	400
LINEARITY LIMIT	8.0
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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ALT/TGP UV

TEMPERATURA: 37°C

TESTNAME	ALT/GPT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	1000
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	-1746

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TEMPERATURA 37°C

TESTNAME	AMICNPG
DECIMAL	0
MEASURE UNIT	U/l
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0.0
MAXIMUM	500
LINEARITY LIMIT	2000
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	3953

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AST/TGO UV

TEMPERATURA: 37°C

TESTNAME	AST/GOT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	1000
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	-1746

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

TESTNAME	CÁLCIO
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	3
REAGENT VOLUME (µL)	
REAGENT 1	300
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1000
LINEARITY LIMIT	20
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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CK - NAC

TEMPERATURA: 37°C

TESTNAME	CK
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	500
LINEARITY LIMIT	1700
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.40
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	8095

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TEMPERATURA: 37°C

TESTNAME	CLORETOS
DECIMAL	0
MEASURE UNIT	mEq/L
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	
LINEARITY LIMIT	125
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0,30
MIX 2	
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	2
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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COLESTEROL 250

TEMPERATURA: 37°C

TESTNAME	COLT
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.40
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEMPERATURA: 37°C

TESTNAME	COLT
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	800
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0,00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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COLESTEROL HDL

TEMPERATURA: 37°C

TESTNAME	CHDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	10
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	400
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

TESTNAME	CREA
DECIMAL	0.0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	500
LINEARITY LIMIT	10.0
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	15
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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DHL - UV

TEMPERATURA: 37°C

TEST NAME	LDH
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	800
MAXIMUM	2000
LINEARITY LIMIT	1257
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 8109

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TEMPERATURA: 37°C

TESTNAME	FALC
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	690
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	2764

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FOSFATO - UV

TEMPERATURA: 37°C

TESTNAME	FÓSFORO
DECIMAL	0.0
MEASURE UNIT	mg/dL
REACTION TYPE	END
POINT	
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1500
LINEARITY LIMIT	20
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

TESTNAME	GGT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	750
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	1111

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GLICOSE ENZ. LÍQUIDA

TEMPERATURA: 37°C

TESTNAME	GLIC
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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MICROPROTE pirogalol

TEMPERATURA: 37°C

TESTNAME	PROTEINÚRIA
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	6
REAGENT VOLUME (µL)	
REAGENT 1	300
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	400
LINEARITY LIMIT	300
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEMPERATURA: 37°C

TEST NAME	PROT.TOT
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	250
LINEARITY LIMIT	12.0
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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TRIGLICÉRIDES 120

TEMPERATURA: 37°C

TESTNAME	TGR
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	1000
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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TEMPERATURA: 37°C

TESTNAME	TGR
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	800
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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URATO 160

TEMPERATURA: 37°C

TESTNAME	AU
DECIMAL	0.0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	14.0
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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URÉIA - UV

TEMPERATURA: 37°C

TEST NAME	UREUV
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	1000
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	10
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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