
LABMAX FLEX CUBETA HELLMA

LANÇAMENTOS!

**CÁLCIO ARSENAZO
MICROPROTE pirogalol**

ÍNDICE

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ALBUMINA

EDIT METHODS					
					Method n. 22
Type	Endpoint				Albumina
Filter(s)	630 nm				Code. ALB
Sample Volume (µL)	3				
Reagent 1 Volume (µL)	750				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	120				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	@	@			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	8.0		N. of decimais	2	
R. Blank Value min	0.00	max	0.2		
Normal Value min	3.50	max	5.50		g/dL
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

ALT/TGP UV

EDIT METHODS									
									Method n. 9
Type	Kinetic								ALT GPT
Filter(s)	340 nm								Code. ALT
Sample Volume	(μ L)	50							
Reagent 1 Volume	(μ L)	500							
Reagent 2 Volume	(μ L)	0		Reagent 3 Volume	(μ L)	0			
Number of Washes	2		Only final washes						
Stabilization Time	(s)	24		Reading Time	(s)	30	Cicle	60	
R1 Incubation Time	(s)	60							
Position of reagent	#								
Number of standards	#			Number. of controls	#				
Position of standard	1			Position of control	7				
Interpolation	Linear regression								
Curve O.D.									
Val. Std	#		#						
Substrate depletion	1.000								
Dilution Ratio in Rep	.5000			Sample dilution	2				
Linearity	350			N. of decimais	0				
R. Blank Value	min	0.8	max	2.200					
Normal Value	min	8.00	max	42.00				U.I./L	
>Blanks are used in calculations									
>The calibration is KEPT in memory									

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REVISÃO: 02 (07/2009)

AMILASE CNPG

EDIT METHODS							
							Method n. 13
Type	Kinetic						Amilase
Filter(s)	405 nm						Code. AMI
Sample Volume (µL)	10						
Reagent 1 Volume (µL)	500						
Reagent 2 Volume (µL)	0			Reagent 3 Volume (µL)	0		
Number of Washes	2		Only final washes				
Stabilization Time (s)	13		Reading Time (s)	17		Cicle	36
R1 Incubation Time (s)	60						
Position of reagent	#						
Number of standards	#			Number. of controls	#		
Position of standard	1			Position of control	7		
Interpolation	Linear regression						
Curve O.D.							
Val. Std	#		#				
Substrate depletion	1.000						
Dilution Ratio in Rep	.5000			Sample dilution	2		
Linearity	2000			N. of decimais	0		
R. Blank Value min	0.0		max	0.200			
Normal Value min	25.0		max	125.0		U/L	
>Blanks are used in calculations							
>The calibration is KEPT in memory							

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REVISÃO: 02 (07/2009)

AST/TGO UV

EDIT METHODS							
							Method n. 10
Type	Kinetic						AST/ TGO
Filter(s)	340 nm						Code. AST
Sample Volume	(μ L)	50					
Reagent 1 Volume	(μ L)	500					
Reagent 2 Volume	(μ L)	0		Reagent 3 Volume	(μ L)	0	
Number of Washes	2		Only final washes				
Stabilization Time	(s)	24		Reading Time	(s)	30	Cicle 60
R1 Incubation Time	(s)	60					
Position of reagent	#						
Number of standards	#		Number. of controls		#		
Position of standard	1		Position of control		7		
Interpolation	Linear regression						
Curve O.D.							
Val. Std	#	#					
Substrate depletion	1.000						
Dilution Ratio in Rep	.5000		Sample dilution		2		
Linearity	350		N. of decimais		0		
R. Blank Value	min	0.8	max	2.200			
Normal Value	min	12.0	max	42.00		U.I./L	
>Blanks are used in calculations							
>The calibration is KEPT in memory							

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CÁLCIO ARSENAZO

EDIT METHODS

					Method n. 4	
Type	Endpoint				Cálcio	
Filter(s)	630 nm				Code. CAR	
Sample Volume (µL)	5					
Reagent 1 Volume (µL)	500					
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0		
Number of Washes	2		Only final washes			
Stabilization Time (s)	4		Reading Time (s)	3		Cicle 12
R1 Incubation Time (s)	180					
Position of reagent	#					
Number of standards	#		Number. of controls	#		
Position of standard	1		Position of contrai	7		
Interpolation	Linear regression					
Curve O.D.						
Val. Std	#		#			
Dilution Ratio in Rep	.5000		Sample dilution	2		
Linearity	20.0		N. of decimais	2		
R. Blank Value min	0.090		max	0.300		
Normal Value min	4.6		max	11.0		mg/dL
>Blanks are used in calculations						
>The calibration is KEPT in memory						

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REVISÃO: 02 (07/2009)

CK-NAC

EDIT METHODS					
					Method n. 25
Type	Kinetic				CK - NAC
Filter(s)	340 nm				Code. CK
Sample Volume (µL)	10				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0	Reagent 3 Volume (µL)	0		
Number of Washes	2	Only final washes			
Stabilization Time (s)	13	Reading Time (s)	17	Cicle	36
R1 Incubation Time (s)	120				
Position of reagent	#				
Number of standards	#	Number. of controls	#		
Position of standard	1	Position of control	7		
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Substrate depletion	1.5				
Dilution Ratio in Rep	.5000	Sample dilution	2		
Linearity	1700	N. of decimais	0		
R. Blank Value min	0.0	max	0.700		
Normal Value min	15	max	110	U.I./L	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

CLORETOS COLORIMÉTRICO

EDIT METHODS					
					Method n. 7
Type	Endpoint				Cloretos
Filter(s)	505nm				Code. CLO
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	120				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of contrai	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#		#		
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	125.0		N. of decimais	2	
R. Blank Value min	0.0		max	0.050	
Normal Value min	96		max	105	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

COLESTEROL 250

EDIT METHODS							
							Method n. 16
Type	Endpoint						Colesterol
Filter(s)	505 nm						Code. COE
Sample Volume (µL)	5						
Reagent 1 Volume (µL)	500						
Reagent 2 Volume (µL)	0			Reagent 3 Volume (µL)	0		
Number of Washes	2		Only final washes				
Stabilization Time (s)	4		Reading Time (s)	3		Cicle	12
R1 Incubation Time (s)	600						
Position of reagent	#						
Number of standards	#			Number. of controls	#		
Position of standard	1			Position of contrai	7		
Interpolation	Linear regression						
Curve O.D.							
Val. Std	#		#				
Dilution Ratio in Rep	.5000			Sample dilution	2		
Linearity	500			N. of decimais	0		
R. Blank Value	min	0.00		max	0.200		
Normal Value	min	50		max	200		mg/dL
>Blanks are used in calculations							
>The calibration is KEPT in memory							

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REVISÃO: 02 (07/2009)

COLESTEROL ENZ. LÍQUIDO**EDIT METHODS**

					Method n. 17	
Type	Endpoint				Colesterol	
Filter(s)	505 nm				Code. COL	
Sample Volume	µL	5				
Reagent 1 Volume	µL	500				
Reagent 2 Volume	µL	0		Reagent 3 Volume	µL	0
Number of Washes	2		Only final washes			
Stabilization Time	(s)	4		Reading Time	(s)	3
R1 Incubation Time	(s)	300		Cicle	12	
Position of reagent	#					
Number of standards	#		Number. of controls	#		
Position of standard	1		Position of contrai	7		
Interpolation	Linear regression					
Curve O.D.						
Val. Std	#	#				
Dilution Ratio in Rep	.5000		Sample dilution	2		
Linearity	800		N. of decimais	0		
R. Blank Value	min	0.00		max	0.200	
Normal Value	min	50		max	200	
>Blanks are used in calculations						
>The calibration is KEPT in memory						

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REVISÃO: 02 (07/2009)

COLESTEROL HDL

EDIT METHODS					
					Method n. 18
Type	Endpoint				Col. HDL
Filter(s)	505 nm				Code. HDL
Sample Volume (µL)	15				
Reagent 1 Volume (µL)	600				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	300				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	@	@			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	200		N. of decimais	0	
R. Blank Value min	0.00	max	0.200		
Normal Value min	40	max	60	mg/dL	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

CREATININA

<u>EDIT METHODS</u>						
						Method n. 23
Type	Fixed Time					Creatinina
Filter(s)	505 nm					Code. CRE
Sample Volume (µL)	50					
Reagent 1 Volume (µL)	500					
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0		
Number of Washes	3		Only final washes			
Stabilization Time (s)	13		Reading Time (s)	60		Cicle 84
R1 Incubation Time (s)	12					
Position of reagent	#					
Number of standards	#		Number. of controls	#		
Position of standard	1		Position of control	7		
Interpolation	Linear regression					
Curve O.D.						
Val. Std	# #					
Substrate depletion	1.2					
Dilution Ratio in Rep	.5000		Sample dilution	2		
Linearity	10.0		N. of decimais	2		
R. Blank Value min	0.00		max	0.50		
Normal Value min	0.5		max	1.2 mg/dL		
>Blanks are used in calculations						
>The calibration is KEPT in memory						

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REVISÃO: 02 (07/2009)

DHL- UV

EDIT METHODS							
							Method n. 27
Type	Kinetic						DHL UV
Filter(s)	340 nm						Code. DHL
Sample Volume (µL)	10						
Reagent 1 Volume (µL)	500						
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0			
Number of Washes	2		Only final washes				
Stabilization Time (s)	13		Reading Time (s)	17		Cicl	36
R1 Incubation Time (s)	60						
Position of reagent	#						
Number of standards	#		Number. of controls	#			
Position of standard	1		Position of control	7			
Interpolation	Linear regression						
Curve O.D.							
Val. Std	#		#				
Substrate depletion	1.000						
Dilution Ratio in Rep	.5000		Sample dilution	2			
Linearity	1257		N. of decimais	0			
R. Blank Value min	0.0		max	0.700			
Normal Value min	207		max	414		U.I./L	
>Blanks are used in calculations							
>The calibration is KEPT in memory							

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FOSFATASE ALCALINA CINÉTICA

EDIT METHODS

					Method n. 12	
Type	Kinetic				Fosf. Alcalina	
Filter(s)	405 nm				Code. FAL	
Sample Volume (µL)	10					
Reagent 1 Volume (µL)	500					
Reagent 2 Volume (µL)			Reagent 3 Volume (µL)	0		
Number of Washes	2		Only final washes			
Stabilization Time (s)	157		Reading Time (s)	17		Cicle 180
R1 Incubation Time (s)	24					
Position of reagent	#					
Number of standards	#		Number. of controls	#		
Position of standard	1		Position of control	7		
Interpolation	Linear regression					
Curve O.D.						
Val. Std	#		#			
Substrate depletion	1.7					
Dilution Ratio in Rep	.5000		Sample dilution	2		
Linearity	690		N. of decimais	0		
R. Blank Value min	0.0		max	0.700		
Normal Value min	100		max	290		U.I./L
>Blanks are used in calculations						
>The calibration is KEPT in memory						

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REVISÃO: 02 (07/2009)

FOSFATO-UV

EDIT METHODS							
							Method n. 2
Type	Endpoint						Fosfato
Filter(s)	340 nm						Code. FOS
Sample Volume (µL)	5						
Reagent 1 Volume (µL)	500						
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0			
Number of Washes	2		Only final washes				
Stabilization Time (s)	4		Reading Time (s)	3		Cicle	12
R1 Incubation Time (s)	300						
Position of reagent	#						
Number of standards	#		Number. of controls	#			
Position of standard	1		Position of control	7			
Interpolation	Linear regression						
Curve O.D.							
Val. Std	#		#				
Dilution Ratio in Rep	.5000		Sample dilution	2			
Linearity	16.0		N. of decimais	2			
R. Blank Value min	0.4		max	0.9			
Normal Value min	2.5		max	4.5			
>Blanks are used in calculations							
>The calibration is KEPT in memory							

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REVISÃO: 02 (07/2009)

GAMA GT CINÉTICA

EDIT METHODS					
					Method n. 11
Type	Kinetic				Gama GT
Filter(s)	405 nm				Code. GGT
Sample Volume (µL)	50				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0	Reagent 3 Volume (µL)	0		
Number of Washes	2	Only final washes			
Stabilization Time (s)	85	Reading Time (s)	17	Cicl	108
R1 Incubation Time (s)	60				
Position of reagent	#				
Number of standards	#	Number. of controls	#		
Position of standard	1	Position of control	7		
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Substrate depletion	1.3				
Dilution Ratio in Rep	.5000	Sample dilution	2		
Linearity	300	N. of decimais	1		
R. Blank Value	min 0.4	max 1.1			
Normal Value	min 10	max 60	U.I./L		
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

GLICOSE ENZ. LÍQUIDA

<u>EDIT METHODS</u>					
					Method n. 28
Type	Endpoint				Glicose
Filter(s)	505 nm				Code. GEL
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	300				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	600		N. of decimais	0	
R. Blank Value min	0.00	max	0.150		
Normal Value min	70	max	99 mg/dL		
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

GLUCOX 500

EDIT METHODS					
					Method n. 29
Type	Endpoint				Glicose
Filter(s)	505 nm				Code. GLU
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)			Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	600				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	600		N. of decimais	0	
R. Blank Value min	0.00	max	0.		
Normal Value min	70	max	99		mg/dL
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

MICROPROTE pirogalol

EDIT METHODS					
					Method n. 30
Type	Endpoint				Microprote
Filter(s)	600 nm				Code. MIC
Sample Volume (µL)	10				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	600				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of contrai	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	100	#			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	300		N. of decimais	2	
R. Blank Value min	0.00	max	0.200		
Normal Value min	7	max	150		mg/dL
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

PROTEÍNAS TOTAIS

<u>EDIT METHODS</u>					
					Method n. 31
Type	Endpoint				Proteínas
Filter(s)	505 nm				Code. PRO
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	300				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#		#		
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	12		N. of decimais	0	
R. Blank Value min	0.00		max	0.900	
Normal Value min	6.0		max	8.0 g/dL	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

TRIGLICÉRIDES 120

<u>EDIT METHODS</u>					
					Method n. 15
Type	Endpoint				Trigl. 120
Filter(s)	505nm				Code. TRI
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	600				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#		#		
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	1000		N. of decimais	0	
R. Blank Value	min	0.00	max	0.500	
Normal Value	min	1	max	200 mg/dL	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

TRIGLICÉRIDES ENZ. LÍQUIDO

<u>EDIT METHODS</u>					
					Method n. 16
Type	Endpoint				Trigl. Enz. Líq.
Filter(s)	505nm				Code. TRE
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	300				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	800		N. of decimais	0	
R. Blank Value	min	0.00	max	0.600	
Normal Value	min	1	max	200 mg/dL	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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REVISÃO: 02 (07/2009)

URATO 160

<u>EDIT METHODS</u>					
					Method n. 21
Type	Endpoint				Urato 160
Filter(s)	505 nm				Code. URA
Sample Volume (µL)	15				
Reagent 1 Volume (µL)	600				
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0	
Number of Washes	2		Only final washes		
Stabilization Time (s)	4		Reading Time (s)	3	Cicle 12
R1 Incubation Time (s)	600				
Position of reagent	#				
Number of standards	#		Number. of controls	#	
Position of standard	1		Position of control	7	
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#		#		
Dilution Ratio in Rep	.5000		Sample dilution	2	
Linearity	14		N. of decimais	2	
R. Blank Value min	0.00		max	0.3	
Normal Value min	1.6		max	7.0 mg/dL	
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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URATO ENZ. LÍQUIDO

EDIT METHODS

					Method n. 33	
Type	Endpoint				Urato Enz. Líquido	
Filter(s)	505 nm				Code. URL	
Sample Volume (µL)	15					
Reagent 1 Volume (µL)	600					
Reagent 2 Volume (µL)	0		Reagent 3 Volume (µL)	0		
Number of Washes	2		Only final washes			
Stabilization Time (s)	4		Reading Time (s)	3		Cicle 12
R1 Incubation Time (s)	600					
Position of reagent	#					
Number of standards	#		Number. of controls	#		
Position of standard	1		Position of control	7		
Interpolation	Linear regression					
Curve O.D.						
Val. Std	#		#			
Dilution Ratio in Rep	.5000		Sample dilution	2		
Linearity	14		N. of decimals	2		
R. Blank Value min	0.00		max	0.3		
Normal Value min	1.6		max	7.0		mg/dL
>Blanks are used in calculations						
>The calibration is KEPT in memory						

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REVISÃO: 02 (07/2009)

UREIA-UV

EDIT METHODS					
					Method n. 32
Type	Fixed Time				Uréia
Filter(s)	340 nm				Code. URE
Sample Volume (µL)	5				
Reagent 1 Volume (µL)	500				
Reagent 2 Volume (µL)	0	Reagent 3 Volume (µL)	0		
Number of Washes	2	Only final washes			
Stabilization Time (s)	13	Reading Time (s)	17	Cicle	36
R1 Incubation Time (s)	60				
Position of reagent	#				
Number of standards	#	Number. of controls	#		
Position of standard	1	Position of control	7		
Interpolation	Linear regression				
Curve O.D.					
Val. Std	#	#			
Substrate depletion	1.2				
Dilution Ratio in Rep	.5000	Sample dilution	1		
Linearity	200	N. of decimais	1		
R. Blank Value	min 0.4	max 0.9			
Normal Value	min 2.5	max 4.5	mg/dL		
>Blanks are used in calculations					
>The calibration is KEPT in memory					

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